

Beskrivelse	Insertion Tube ø (mm)	Working Channel ø (mm)	Working Length (mm)	View Angle (°)
<b>Gastroskoper</b>				
EG27-i10 Video Gastroscope i10 Series	9	2,8	1050	140
EG29-i10 Video Gastroscope i10 Series	9,8	3,2	1050	140
EG34-i10 Video Gastroscope	11,6	3,8	1050	140
EG-2990Zi Video Gastroscope 90i Series MagniView	9,8	2,8	1050	140
EG-2990K Video Gastroscope 90K Series	9,8	2,8	1050	140
EG-2490K Video Gastroscope 90K Series	8	2,4	1050	140
EG-2790K Video Gastroscope 90K Series	9	2,8	1050	140
EG-2990i Video Gastroscope 90i Series	9,8	2,8	1050	140
EG-3890TK Video Gastroscope 90K Series With Double Channel	12,8	3,8/2,8	1050	140
NEW! EG29-i10c IMAGINA Video Gastroscope i10c Series	9,8	3,2	1050	140