

Beskrivelse	Insertion Tube \varnothing (mm)	Working Channel \varnothing (mm)	Working Length (mm)	View Angle ($^{\circ}$)
Bronkoskoper				
EB15-J10 Video Bronchoscope J10 Series	5,2	2	600	120
EB19-J10 Video Bronchoscope J10 Series	6,4	2,8	600	120
EB-1170K Video Bronchoscope 70K Series	3,7	1,2	600	120
EB-1970TK Video Bronchoscope 70TK Series	6,3	3,2	600	120
EB-1990i Video Bronchoscope 90i Series	6,4	1,2	600	120
FB-18V Fiber Bronchoscope V Series	6	2,7	600	120
FB-19TV Fiber Bronchoscope V Series	6,2	3,1	600	120
FB-18RBS Fiber Bedside Bronchoscope	6	2,55	600	100
FB-15RBS Fiber Bedside Bronchoscope	4,9	1,95	600	100
FB-15V Fiber Bronchoscope V Series	4,9	2,1	600	120
FB-10V Fiber Bronchoscope V Series	3,5	1,2	600	120
FB-8V Fiber Bronchoscope V Series	2,8	1,2	600	100