

Produktgruppe	Beskrivelse	Insertion Tube ø (mm)	Working Channel ø (mm)	Working Length (mm)	View Angle (°)
Ultralydskoper					
Gastroskop	EG34-J10U Video Ultrasound Gastroscope J10 Series	11,6	2,8	1250	120
Gastroskop	EG36-J10UR Video Ultrasound Gastroscope J10 Series	12,1	2,4	1250	140
Gastroskop	EG38-J10UT Video Ultrasound Gastroscope J10 Series	12,8	4,0	1250	120
Gastroskop	EG-3270UK Video Ultrasound Gastroscope 70K Series	10,8	2,8	1250	120
Gastroskop	EG-3670URK Video Ultrasound Gastroscope 70K Series	12,1	2,4	1250	140
Gastroskop	EG-3870UTK Video Ultrasound Gastroscope 70K Series	12,8	3,8	1250	120
Bronkoskop	EB19-J10U Video Ultrasound Bronchoscope J10 Series	6,3	2,2	600	100
Bronkoskop	EB-1970UK Video Ultrasound Bronchoscope 70K Series	6,3	2	600	120