

Beskrivelse	Insertion Tube \varnothing (mm)	Working Channel \varnothing (mm)	Working Length (mm)	View Angle ($^{\circ}$)
Naso-Pharyngo-Laryngoskoper (\varnothingNH)				
VNL9-CP VIVIDEO Video Naso-Pharyngo- Laryngoscope CP Series	3,3	0	300	90
VNL8-J10 HD Video J10 Series Naso-Pharyngo-Laryngoscope	2,9	0	300	80
VNL11-J10 HD Video J10 Series Naso-Pharyngo-Laryngoscope	3,6	0	300	80
VNL15-J10 HD Video J10 Series Naso-Pharyngo-Laryngoscope	4,9	2	300	80
VNL-1190STK Video 90K Series Naso-Pharyngo-Laryngoscope	3,7	0	300	80
VNL-1590STi HD Video Naso-Pharyngo-Laryngoscope	5,1	0	300	100